FLUORION

nCoV-19 QLP 2.1

Real-Time PCR Kiti

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Transfusion Transmitted Virus (TTV)	-	+
Human Papilloma Virus (HPV)	-	+
Eppstein Barr Virus (EBV)	-	+
Hepatitis B Virus (HBV)	-	+
Herpes Simplex Virus 1 (HSV 1)	-	+
Herpes Simplex Virus 2 (HSV 2)	-	+
Influenza A Virus	-	+
Influenza B Virus	-	+
Hepatitis D Virus	-	+
Rhinovirus	-	+
Cytomegalo Virus	-	+
Respiratory Syncytial Virus	-	+
Mycobacterium Tuberculosis	-	+
H1N1 Virus	-	+

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13. REFERANSLAR

Espy MJ1, Uhl JR, Sloan LM, Buckwalter SP, Jones MF, Vetter EA, Yao JD, Wengenack NL, Rosenblatt JE, Cockerill FR 3rd, Smith TF., (2006), Real-time PCR in clinical microbiology: applications for routine laboratory testing. *Clin Microbiol Rev.* 19(1),165-256.

WHO. Infection prevention and control during health care when novel coronavirus disease (COVID-19) is suspected or confirmedexternal icon.29 June 2020.

WHO. Transmission of SARS-CoV-2: implications for infection prevention precautionsexternal icon. 9 July 2020.

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TIONTEK

Table of Contents

	<u>Page</u>
1. Kit Contents	1
2. Storage	1
3. Additional Required Materials and Devices	1
4. Important Notes and Safety Instructions	2
5. Pathogen Information	3
6. Principles of the Procedure	3
7. Product Description	4
8. Intended Use	4
9. Product Use Limitations	5
10. Test Procedure	5
10.1. Sample Collection, Storage and Transport	5
10.2. Interfering Substances	5
10.3. RNA Isolation	6
10.4. Specifications of the Kit Elements	6
10.4.1. Negative Control	6
10.4.2. Master Mix	6
10.4.3. Positive Control	7
10.5. Preparing the PCR	7
10.6. Thermal Protocol	8
10.7. Data Analysis	8
11. Troubleshooting	9
12. Product Specifications	10
12.1. Sensitivity	10
12.2. Cross-Reactivity	10

12.3. Reproducibility	11
13. References	12
14. Ordering Information	12
15. Symbols	13
16. Contact Information	14

1. KIT CONTENTS

FLUORION nCoV-19 QLP 2.1 Real-Time PCR Kit is composed of PCR reagents (Table 1):

Table 1. Kit Contents

TUBE NO.	REAGENT	1000 Rxn	500 Rxn	100 Rxn	25 Rxn
1	dH ₂ O*	(1000 µl)	(1000 µl)	(500 µl)	(500 µl)
2a	Master Mix	(8500 µl)	(4250 µl)	(850 µl)	(212.5 µl)
6a	Positive Control	(600 µl)	(450 µl)	(150 µl)	(150 µl)

*dH₂O is used as Negative Control.

2. STORAGE

All of the **FLUORION** nCoV1-9 QLP 2.1 Real-Time PCR Kit components should be stored at -20°C or below. Repeated thawing and freezing (>2x) should be avoided since it may reduce sensitivity of the kit. If the components are to be used in small amounts, they should be frozen in aliquots.

While preparing the PCR the components should not be exposed to room temperature for more than 10-15 min. and the detection mix components should not be exposed to light (more than 1-2 min).

The components maintain their stability until the expiry dates on the labels if they are stored at proper conditions.

3. ADDITIONAL REQUIRED MATERIALS AND DEVICES

- Fluorion Detection System (Iontek), Rotorgene (Qiagene), Tianlong, LightCycler 480 (Roche), Biometrics LIGHTGene 40G/60G, Anitoo Maverick, Heal Force x960 Real-Time PCR System, Quant Studio 5 (Thermo Fisher), Bio-Rad CFX96
- PCR strips or tubes (0,2 ml, DNase, Rnase-free polypropylene tubes or strips with flat cap).
- Fluorion i12, i24 Automated Isolation System / i12 Viral RNA Extraction Kit (Iontek), ZipPrime Viral Nucleic Acid Isolation Kit, Qiagen MinElute Virus Spin Kit, Smart NAT Viral Transfer Solution
- Deep freezer (-20°C or below)
- Heat Block
- Vortex mixer
- Desktop centrifuge with rotor for 2 ml. microcentrifuge tubes

- Adjustable micropipettes
- Sterile micropipette tips with filters
- Sterile 1.5 or 2 ml. microcentrifuge tubes
- Disposable gloves (Powder-free)

4. IMPORTANT NOTES AND SAFETY INSTRUCTIONS

Caution!:

- The product should be delivered on dry ice. Check for presence of dry ice upon arrival, and do not use products that have been delivered without it.
- Check for the expiry dates on the box and tube labels, upon arrival. Do not use expired products or components.

Attention should be paid also to the following points:

- Sterile micropipette tips with filters, and sterile microcentrifuge tubes should be used
- Before starting a test procedure all components should be thoroughly thawed at room temperature. After thawing, all components should be mixed and centrifuged briefly, to ensure homogeneity prior to use.
- The kit components should be kept on ice or a cooling block until the reaction is prepared, and they should be immediately returned to -20°C or below.
- PCR and nucleic acid isolation should be performed in different compartments. Samples should be stored separately to avoid contact with the kit components.

The safety instructions below should be followed:

- Samples (including the positive control) should be handled with extreme caution: Physical contact with pathogens should be avoided by; wearing lab coats and gloves, no allowance for eating or drinking within the workspace, prevention of unauthorized individuals' access to the working area.
- Pathogenic wastes (biological samples and material contaminated with them) produced during the nucleic acid isolation step, should be discarded into medical waste and disposed safely.

- In case of a spill the contaminated area must be cleaned with %5 hypochloride solution and exposed to UV for 15 minutes to remove hazardous nucleic acids or pathogens.
- Pathogen information should be reviewed for the health risks.

5. PATHOGEN INFORMATION

COVID-19, which belongs to coronaviruses family, was first identified in Wuhan, China in December 2019. COVID-19 is caused by the virus severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2), a highly spreadable new virus in humans causing respiratory illness. Transmission from person-to-person happens through respiratory droplets which are released when someone with COVID-19 sneezes, coughs, or talks. Moreover, recent data suggest that people who don't show any or mild symptoms can also transmit COVID-19 through those droplets. Covid-19 symptoms include fever, cough, shortness of breath, headache, fatigue, sore throat [2,3].

6. PRINCIPLES OF THE PROCEDURE

FLUORION nCoV-19 QLP 2.1 Real-Time PCR Kit is based on the real-time RT-PCR principle. The pathogen is detected using fluorescent dyes that are incorporated into oligonucleotide probes. The assay utilizes the 5' exonuclease activity of Taq Polymerase to cleave a dual-labeled fluorescent hybridization probe during the extension phase of PCR. The probe is labeled at the 5' end with a fluorescent 'reporter' molecule, and at the 3' end with another fluorescent molecule that acts as a 'quencher' for the 'reporter'. When the two fluorophores are in close proximity and the reporter is excited by light, no reporter fluorescence can be detected. During the elongation step of PCR, Taq Polymerase encounters and cleaves the probe bound to the template. As the reporter is freed from the suppressing effect of the quencher, fluorescence signal can be detected. The fluorescence generated by the reporter increases as the PCR product is accumulated; the point at which the signal rises above background level and becomes distinguishable, is called the threshold cycle (C_T). There is a linear relationship between the log of the starting amount of a template and its threshold cycle, thus starting amount of unknown templates can be determined using standard curves constructed using C_T values of the known starting amounts of target templates.

During thermal cycling, reverse transcription and PCR takes place in the same tube. All reagents required for both reactions are added at the beginning, so there is

no need to open the tube after the reverse transcription reaction has been finished. The Taq Polymerase inside the reaction mix is a hot start enzyme which is provided in an inactive state. It has no enzymatic activity at ambient temperatures. The enzyme remains completely inactive during the reverse transcription reaction and does not interfere with it. This prevents the formation of misprimed products and primer-dimers during reaction setup, reverse transcription, and the first denaturation step. The enzyme is activated by a 2-minute, 95 °C incubation step, which is integrated in the thermal cycling program. The hot start also inactivates the reverse transcriptase enzymes, ensuring temporal separation of reverse transcription and PCR, and allowing both steps to be performed sequentially in a single tube.

In contrast with conventional PCR, Real-Time PCR obviates the need for further analysis methods like gel electrophoresis, whereby minimizing the risk of contamination.

7. PRODUCT DESCRIPTION

FLUORION nCoV-19 QLP 2.1 Real-Time PCR Kit detects SARS-CoV-2 RNA in oropharyngeal and nasopharyngeal swab samples. The analytic sensitivity is 10 copies/reaction. N1 and N2 genes are amplified using sequence-specific primers and the detection is done using the FAM filter pair. To control PCR inhibition an endogenous internal control is incorporated into the system. Endogenous internal control is detected with the HEX filter pair during data collection.

This kit has been validated with Fluorion Detection Systems (Iontek), Rotorgene (Qiagen), LightCycler 480 (Roche), Gentier 48E (Tianlong), Biometrics LIGHTGene 40G/60G, Anitox Maverick, Heal Force x960 Real-Time PCR System, Quant Studio 5 (Thermo Fisher) Real-Time PCR Systems and i-Series extraction systems, ZipPrime Viral Nucleic Acid Isolation Kit, Qiagen MinElute Virus Spin Kit, Smart NAT Viral Transport Solution.

8. INTENDED USE

The **FLUORION** nCoV-19 QLP 2.1 Real Time PCR Kit is an in vitro diagnostic kit designed for the quantitative measurement of SARS-COV-2 RNA in swab samples. The results from the **FLUORION** nCoV-9 QLP 2.1 Real Time PCR Kit must be interpreted within the context of all relevant clinical and laboratory findings.

The **FLUORION** nCoV-19 QLP 2.1 Real Time PCR Kit is not intended for use as a screening test for the presence of SARS-COV-2 RNA in blood or blood products.

9. PRODUCT USE LIMITATIONS

- All the components may exclusively be used for in vitro diagnostics.
- This product should be used in accordance with this user manual.
- This product is to be used by personnel specially trained to perform in vitro diagnostic procedures and applying molecular biology technics.
- Rare mutations within the regions of the viral genome covered by the primers or probe in **FLUORION** nCoV-19 QLP 2.1 Real Time PCR Kit may result in failure detection of virus. **FLUORION** nCoV-19 QLP 2.1 Real Time PCR Kit primers and probe are designed to target the most conserved region of the N1 and N2 genes. **FLUORION** nCoV-19 QLP 2.1 Real Time PCR Kit detects the SARS-CoV-2 virus independently of these variants, as the mutations of the B.1.1.7 (UK), B.1.351 (South Africa), P.1 (Brazil), B.1.525 (New York), B.1.429 (California) and similar variants do not overlap with the primer / probe region.
- During the early phase of infection (window period) **FLUORION** nCoV-19 QLP 2.1 Real Time PCR Kit may result in failure of detection of the virus.

10. TEST PROCEDURE

10.1. Sample Collection, Storage and Transport

Nasopharyngeal and oropharyngeal swab samples must be collected according to the relevant sampling protocols.

The samples should be transported frozen in dry ice in containers with capacity to resist pressure. Transportation should be done according to local and national regulations for pathogen material transport.

10.2. Interfering Substances

To avoid possible influences on PCR:

- Samples which have been collected in tubes that contain anticoagulants such as heparin,
- Hemolytic samples,
- Samples of heparinised patients,
- Samples of patients with elevated levels of bile salts, bilirubin or lipids

must not be used.

10.3. RNA Isolation

FLUORION nCoV-19 QLP 2.1 Real-Time PCR Kit is optimized with **FLUORION** i12, i24 (Iontek), ZipPrime Viral Nucleic Acid Isolation Kit, Qiagen MinElute Virus Spin Kit automated and manual extraction systems and corresponding kits. The RNA isolation should be performed according to the manufacturer's instructions.

10.4. Specifications of the Kit Elements

10.4.1. Negative Control

PCR-grade dH₂O included in the kit is used as negative control.

10.4.2. Master Mix

HotStarTaq DNA Polymerase: HotStarTaq DNA Polymerase is a modified form of a recombinant 94 kDa DNA polymerase, originally isolated from *Thermus aquaticus*, cloned into *E. coli*. The enzyme is provided in an inactive form. It is activated by a 2-minute 95 °C incubation step. This prevents the formation of misprimed products and primer-dimers during reaction setup and the first denaturation step, leading to high PCR specificity and accurate quantification.

Probe RT-PCR Buffer: Contains Tris-Cl, KCl, (NH₄)₂SO₄, 8 mM MgCl₂, pH 8.7 (20°C).

dNTP Mix: Contains ultrapure quality dATP, dGTP, dCTP ve dTTP/dUTP.

RT Mix: Contains an optimized mixture of Reverse Transcriptases, expressed in *E. coli*. The optimal working temperature of the enzymes is 46 °C.

Detection Mix combination: Contains SARS-COV-2 -specific forward and reverse primers and a dual-labeled probe (FAM-BHQ1), RNaseP specific dual-labeled probe (HEX-BHQ2).

An endogenous internal control is included in the kit to control isolation and PCR inhibition. The internal control should have a Ct value between 20 < Ct < 40. A larger Ct value or a smaller RFU indicates a problem in isolation or PCR inhibition. In this case the isolation and PCR should be repeated. In samples containing a high viral load, the internal control can be suppressed, and no increase of the signal is detected. Please use the table (Table 2) below for the interpretation of internal control data:

Table 2. Data Interpretation

SARS-COV-2 (FAM)	Internal Control (HEX)	Interpretation
+	+	Sample positive
+	-	Sample positive
-	+	Sample negative
-	-	Repeat the test!

10.4.3. Positive Control

The positive control contains SARS-COV-2 and RNaseP nucleic content. It can be included in the PCR to test the efficiency of the PCR exclusively. The threshold cycle for the positive control is given in the acceptance criteria table (Table 9). Threshold cycles higher than the acceptance criteria may indicate an efficiency loss in the reaction.

10.5. Preparing the PCR

Samples, positive control and negative control (PCR-grade water) should be added to the reaction following the amounts listed in Table 3. Please make sure that all the kit components are thawed before use.

Table 3. PCR reaction volumes	
PCR Mix	8.0 µl
Sample RNA Standard Negative/Positive Control	2.0 µl
Total Volume	10.0 µl

Pipette 8 µl of the master mix into the PCR tubes and add 2 µl of RNA (sample/positive control/negative control). Close the caps of the tubes. Make sure that the cap is firmly closed without leaving any opening on the tubes and the solution in each tube is at the bottom of the tube. Centrifuge if necessary.

10.6. Thermal Protocol

The thermal protocol for **FLUORION** nCoV-19 QLP 2.1 Real-Time PCR Kit is composed of a reverse transcription, an initial denaturation for activation the HotStarTaq DNA Polymerase, a two-step amplification cycle and a terminal hold (Table 4). The Real-Time data is collected at the second step of the amplification cycle.

Table 4. PCR Programme

Step	Temp.	Time	Cycles	
1	Reverse Transcription	46°C	300 seconds (5 min.)	1
2	Initial Denaturation	95°C	120 seconds (2 min.)	1
3	Denaturation	98°C	5 seconds	40
	Annealing and Elongation (Data Collection)	58°C	10 seconds	

***Recommended ramp rate is > 4 °C/sec.**

10.7. Data Analysis

Analysis of the results should be performed by trained personnel who received the required training for analysing Real-Time data. We recommend that the results of the tests must be evaluated by an expert clinician, taking the patient’s clinical findings and the results of other tests into consideration.

Table 5 and Table 6 below display the acceptance criteria and data interpretation for **FLUORION** nCoV-19 QLP 2.1 Real-Time PCR Kit.

Table 5. Kit Acceptance Criteria

Component/Parameter	Cycle Threshold (C _T)
Positive Control	20 < Ct < 30
Internal Control	20 < Ct < 40

The following table (Table 6) shows the possible results and their interpretation:

Table 6. Data Interpretation

Signal detected in FAM filter pair	The sample contains SARS-COV-2 RNA, the result is positive	No need to check the internal control since the sample is positive (high positive samples may suppress the signal from the internal control)
No signal in FAM, signal in HEX	The SARS-COV-2 RNA in the sample is not detectable	Signal from HEX filter pair rules out the possibility of PCR inhibition
No signal in FAM and HEX	The diagnosis is inconclusive	No signal in HEX points out to PCR inhibition or to a problem in RNA isolation (See 11. Troubleshooting)

11. TROUBLESHOOTING

Please contact the manufacturer in case of a problem during a run.

Late or no signal from the FAM filter

Wrong thermal protocol is chosen	Make sure that the correct thermal protocol is chosen.
Polymorphism in the primer or probe binding regions	This possibility is very low. If there is an incompatibility with other clinical test results, it is recommended to confirm the results with alternative tests.

No signal from the internal control

Deterioration of the internal control or detection mix 2	Follow the instructions for the storage of kit components (See 2. Storage).
PCR inhibition	Make sure that you use the recommended RNA isolation method (See 10.3 RNA isolation).
RNA lost during isolation	Make sure that you use the recommended RNA isolation method (See 10.3 RNA isolation).

Signal from FAM Filter in the Negative Control

Contamination	Use filter-tips. Repeat PCR with new kit components.
---------------	--

The Threshold is Above Low Signals

The threshold should be adjusted	Pull the threshold down until it cuts the low signals. Avoid the background and the signal from negative control.
----------------------------------	---

Uneven Traces

Baseline cycles should be readjusted	Assing new baseline cycles manually and recalculate threshold cycles.
--------------------------------------	---

12. PRODUCT SPECIFICATIONS

12.1. Sensitivity

Analytical sensitivity may be expressed as the limit of detection: i.e. the smallest amount of the target marker that can be precisely detected. The detection limit of an individual analytical procedure is the lowest amount of nucleic acid in a sample which can be detected but not necessarily quantitated as an exact value. The analytical sensitivity or detection limit for NAT assays is expressed by the 95% positive cut-off value. This is the analyte concentration where 95% of test runs give positive results following serial dilutions of an international reference material.

The analytical detection limit was determined using serial dilutions of the isolated serum standard with the lowest concentration. Each dilution was tested at least 24-times and the results were analyzed by probit method.

The analytical detection limit for **FLUORION** nCoV-19 QLP 2.1 Real-Time PCR Kit was found to be 10 copies/ml ($p=0.05$), which means that there is a 95% probability that 10 copies/ml will be detected.

12.2 Cross-Reactivity

To eliminate potential cross-reactivity, assay design evidence was used. Primer and probe sequences were checked for possible homology to other known pathogen sequences by sequence comparison analysis using database alignment.

Additionally, the below listed pathogens were tested against cross reactivity and none was detected.

Table 7. Cross Reactivity Data

Pathogen	SARS-COV-2 FAM	Internal Control HEX
Parvovirus B19	-	+
Transfusion Transmitted Virus (TTV)	-	+
Human Papilloma Virus (HPV)	-	+
Eppstein Barr Virus (EBV)	-	+
Hepatitis B Virus (HBV)	-	+
Herpes Simplex Virus 1 (HSV 1)	-	+
Herpes Simplex Virus 2 (HSV 2)	-	+
Influenza A Virus	-	+
Influenza B Virus	-	+
Hepatitis D Virus	-	+
Rhinovirus	-	+
Cytomegalo Virus	-	+
Respiratory Syncntial Virus	-	+
Mycobacterium Tuberculosis	-	+
H1N1 Virus	-	+
Crimean Congo Hemorrhagic Fever Virus	-	+

12.3 Reproducibility

Reproducibility data (on C_T value basis) were obtained by the analysis of the positive control provided by the **FLUORION** nCoV-19 QLP 2.1 Real-Time PCR Kit. Test was performed in at least 4 replicates by 3 different operators, on multiple days, using 3 different lots. The resulting data is given in Table 8.

Table 8. Reproducibility Data

SARS-COV-2 Positive Control	Standard deviation	Variance	Coefficient of variation [%]
Intra-assay Variability N=8	0,060	0,004	0,189
Inter-lot Variability N=3	0,220	0,048	0,692
Inter-operator Variability N=3	0,053	0,003	0,167
Total Inter-assay Variability N=6	0,186	0,035	0,586

13. REFERENCES

Espy MJ1, Uhl JR, Sloan LM, Buckwalter SP, Jones MF, Vetter EA, Yao JD, Wengenack NL, Rosenblatt JE, Cockerill FR 3rd, Smith TF., (2006), Real-time PCR in clinical microbiology: applications for routine laboratory testing. *Clin Microbiol Rev.* 19(1),165-256.

WHO. Infection prevention and control during health care when novel coronavirus disease (COVID-19) is suspected or confirmed external icon. 29 June 2020.

WHO. Transmission of SARS-CoV-2: implications for infection prevention precaution external icon. 9 July 2020.

14. ORDERING INFORMATION

Kit Name	Catalog No.	Test Number
Fluorion nCoV-19 QLP 2.1 Real-Time PCR Kit	M1350102-2	50
Fluorion nCoV-19 QLP 2.1 Real-Time PCR Kit	M1350102-3	100

15. SYMBOLS

Expiration Date

Lot/Batch

Catalog number

Temperature limitation

Caution; consult accompanying documents

Manufacturer

In Vitro Diagnostic Medical Device

Consult instructions for use

Contains sufficient for (n) amount tests

16. CONTACT INFORMATION

Sultan Selim Mh. Turan Sk.
No.21/1
34415 Kağıthane İstanbul-TURKEY
Phone: +90 212 481 55 16
Fax: +90 212 481 55 18
e-mail: iontek@iontek.com.tr
www.iontek.com.tr

Trademarks: ®, QuantiTect™, HotStarTaq™, Sensiscript™, Omniscript™, and QIAamp™ are registered trademarks of QIAGEN GmbH, Hilden, Germany. QIAGEN components contained herein are developed and manufactured by QIAGEN GmbH, Hilden, Germany.

Accurun controls and Basematrix Diluent are trademarks of SeraCare Diagnostics.
IONTEK and FLUORION are registered trademarks of Iontek Inc.